

Essure Hysteroscopic Procedure Equipment

MANDATORY

- Hysteroscope**
*Rigid, 5fr working channel, OD 5mm (7mm max),
12-30° lens, continuous flow*
- Sealing Cap**
For working channel
- Light Source and Cable**
- Camera Source and Head**
- Monitor**
- Tall IV Pole**
- Gauze Forceps**
- Bivalve, Open-Side Speculum**
e.g., Collins
- Single Toothed Tenaculum x2**
- Stirrups**
Ski boot style

STAND-BY

- Pressure Bag**
- Cervical Dilation Set**
- Multi-Toothed Graspers**

OPTIONAL

- VCR**

Essure Hysteroscopic Procedure Disposables & Supply Items

MANDATORY

- Essure Systems**
2 per box (plus back-ups)
- DryFlow™ Introducers**
2 per box
- Patient ID Card**
Save for Post-op
- Antibacterial Prep**
- Gauze**
- Warm Saline**
2 or 3 liter bags
- Bifurcated Large Bore Irrigation Tubing**
- Local Anesthetic**
Supplies to administer local anesthetic

OPTIONAL

- Under Buttocks Drape with Fluid Control Pouch**
- Lithotomy Drapes**
- Sterile Bag/Towel**

Conceptus Incorporated
331 East Evelyn Avenue
Mountain View, CA 94041

Tel: 1-877-ESSURE2 or
1-877-377-8732

www.essureMD.com